



**KARMA HEALTHCARE**  
 19255 San Jose Ave.  
 City Of Industry, CA 91748  
**Tel:** 626-581-2235 **Fax:** 626-581-2335  
**Toll Free:** 800-805-2762



QR Reader Smartphone Info



**FLEXX WHEELCHAIR ORDER FORM**

Submitting for  QUOTE  ORDER

**Mark For:** \_\_\_\_\_  
 Date: \_\_\_\_\_  
 Dealer Acct#: \_\_\_\_\_  
 Dealer: \_\_\_\_\_  
 Dealer Contact: \_\_\_\_\_  
 Dealer Address: \_\_\_\_\_  
 Dealer City: \_\_\_\_\_ ST: \_\_\_\_\_ ZIP: \_\_\_\_\_  
 Dealer Phone: (\_\_\_\_) \_\_\_\_\_ FAX: (\_\_\_\_) \_\_\_\_\_  
 Confirmation Email: \_\_\_\_\_

**Purchase Order #:** \_\_\_\_\_  
 Ship To: \_\_\_\_\_  
 Attention: \_\_\_\_\_  
 Address: \_\_\_\_\_  
 Address: \_\_\_\_\_  
 Ship To City: \_\_\_\_\_ ST: \_\_\_\_\_ ZIP: \_\_\_\_\_  
 Ship To Phone: \_\_\_\_\_

*The HCPCS CODES are not intended to be, nor should be considered billing or legal advice. They are provided based on code verification through PDAC or our interpretations of billing. The use of HCPCS CODES does not ensure coverage or payment. For coverage information, please verify the appropriate payer's coverage policy*

**FLEXX Wheelchair:**

*FLEXX Ultra lightweight Adjustable Wheelchair standard features: swing away removable footrest, fixed backrest, quick release 24"x1" solid poly-urethane rear wheels, 8"x1" solid front casters, adjustable front forks, weighs 28.5 lbs, flip back height adjustable armrest, weight capacity 286 lbs.*

BASE FRAME				
✓	SKU#	HCPCS	Description	Price
<input type="checkbox"/>	KM8522Q	K0005	"FLEXX" Wheelchair Diamond Black Color	\$1899.00

  

SEAT WIDTH				
<input type="checkbox"/>	16B		16" wide	Std
<input type="checkbox"/>	18B		18" wide	Std
<input type="checkbox"/>	20B	E2201	20" wide	\$300.00

  

SEAT DEPTH				
<input type="checkbox"/>	D16		16" depth (not available for 20" seat width)	Std
<input type="checkbox"/>	D18		18" depth	Std

  

FRONT SEAT TO FLOOR HEIGHT				
<input type="checkbox"/>	FSH15	K0056	15" seat to floor height	\$200.00
<input type="checkbox"/>	FSH16	K0056	16" seat to floor height	\$200.00
<input type="checkbox"/>	FSH17		17" seat to floor height	Std
<input type="checkbox"/>	FSH18		18" seat to floor height	Std
<input type="checkbox"/>	FSH19		19" seat to floor height	Std

FRONT SEAT TO FLOOR HEIGHT cont.				
✓	SKU#	HCPCS	Description	Price
<input type="checkbox"/>	FSH20		20" seat to floor height	Std
<input type="checkbox"/>	FSH21		21" seat to floor height	Std

REAR SEAT TO FLOOR HEIGHT*Maximum 2" difference from front seat to floor height*				
<input type="checkbox"/>	RSH15	K0056	15" seat to floor height	See FSH
<input type="checkbox"/>	RSH16	K0056	16" seat to floor height	See FSH
<input type="checkbox"/>	RSH17		17" seat to floor height	Std
<input type="checkbox"/>	RSH18		18" seat to floor height	Std
<input type="checkbox"/>	RSH19		19" seat to floor height	Std
<input type="checkbox"/>	RSH20		20" seat to floor height	Std

SEAT UPHOLSTERY				
<input type="checkbox"/>	US-8500		Non-Adjustable	Std
<input type="checkbox"/>	OMITUS-8500		Omit Seat Upholstery	Std

FRONT CASTER SELECTION (pair)				
<input type="checkbox"/>	C51LT		5"x1" composite hub poly urethane flat free	Std
<input type="checkbox"/>	C5EPIC		5"x1" aluminum hub urethane narrow court caster	\$90.00
<input type="checkbox"/>	C61LT		6"x1" composite hub poly urethane flat free	Std
<input type="checkbox"/>	C6EPIC		6"x1" aluminum hub urethane narrow court caster	\$108.00
<input type="checkbox"/>	C6EPICW		6"X1.4" aluminum hub wide soft roll caster	\$132.00
<input type="checkbox"/>	C71LT		7"x1" composite hub poly urethane flat free	Std
<input type="checkbox"/>	C81LT		8"x1" composite hub poly urethane flat free	Std

Caster Size/Front Seat to Floor Selection				
Front Seat to Floor	Caster Size			
	5"	6"	7"	8"
15"	X			
16"	X	X		
17"	X	X		
18"	X	X	X	
19"		X	X	X
20"				X
21"				X

ARM OPTIONS				
<input type="checkbox"/>	NAARMS		No armrest, No arm pad	Std
<input type="checkbox"/>	AR8500SAP		Swing away padded arms with round tube arm pads	N/A
<input type="checkbox"/>	AR8500HASP	E0973	Flip back height adjustable single post arms with ergonomic arm pads	Std

BACK POST				
<input type="checkbox"/>	FBH8500-17		17" w/ergonomic push handles, 8 degree bend	Std
<input type="checkbox"/>	FBH8500-20		20" w/ergonomic push handles, 8 degree bend	Std
<input type="checkbox"/>	ABH8500-1521		15"-21" range, Fixed Back w/ergonomic push handles, 8 degree bend	\$165.00

BACK UPHOLSTERY				
✓	SKU#	HCPCS	Description	Price
<input type="checkbox"/>	UB-8500-FLEXX		Slip on padded back upholstery	Std
<input type="checkbox"/>	OMITUB-8500		Omit back upholstery	Std
<input type="checkbox"/>	ADJTENS8500		Adjustable Tension back upholstery	N/A

CENTER OF GRAVITY				
<input type="checkbox"/>	COGADJ1		Rear wheel center of gravity adjustment: Stable	Std
<input type="checkbox"/>	COGADJ0		Rear wheel center of gravity adjustment: Normal	Std
<input type="checkbox"/>	COGADJ-1		Rear wheel center of gravity adjustment: Active	Std

REAR WHEELS (pair)				
<input type="checkbox"/>	RW-8500Q		24x1 quick release flat free polyurethane wheel assem. with Ergonomic handrim	Std
<input type="checkbox"/>	RW88		24x1 MAG 8 spoke flat free polyurethane wheel assembly	N/A
<input type="checkbox"/>	SPINLX		24x1 Spinergy LX 12 spoke, Pneumatic Kenda Tire and Tube w/ anodized polished handrim	\$895.00
<input type="checkbox"/>	SPINSPOX		24x1 Spinergy SPOX 18 spoke, Pneumatic Kenda Tire and Tube w/anodized polished handrim	\$695.00
<input type="checkbox"/>	SPINWIRE		24x1 Spinergy WIRE 30 spoke, Pneumatic Kenda Tire and Tube w/anodized polished handrim	\$295.00

KENDA TIRE COLOR/SPINERGY SPOKE COLOR (Selection applies only to Spinergy Wheel)				
<input type="checkbox"/>	KENDAYLW		Yellow Kenda tire	Std
<input type="checkbox"/>	KENDARED		Red Kenda tire	Std
<input type="checkbox"/>	KENDABLU		Blue Kenda tire	Std
<input type="checkbox"/>	SPINSPYLW		Yellow Spinergy Spoke color	Std
<input type="checkbox"/>	SPINSPRED		Red Spinergy Spoke color	Std
<input type="checkbox"/>	SPINSPBLU		Blue Spinergy Spoke color	Std
<input type="checkbox"/>	SPINSPBLK		Black Spinergy Spoke color	Std

BACKREST ANGLE				
<input type="checkbox"/>	BAADJ-4		Back angle adjustment -4 degrees	Std
<input type="checkbox"/>	BAADJ0		Back angle adjustment 0 degrees	Std
<input type="checkbox"/>	BAADJ4		Back angle adjustment 4 degrees	Std
<input type="checkbox"/>	BAADJ8		Back angle adjustment 8 degrees	Std
<input type="checkbox"/>	BAADJ12		Back angle adjustment 12 degrees	Std
<input type="checkbox"/>	BAADJ16		Back angle adjustment 16 degrees	Std

HANGAR OPTIONS				
<input type="checkbox"/>	FR8500Q-16		16" 70 degree Swing in/out non tapered footrest (pair)	Std
<input type="checkbox"/>	FR8500Q-18		18" 70 degree Swing in/out non tapered footrest (pair)	Std
<input type="checkbox"/>	FR8500Q-20		20" 70 degree Swing in/out non tapered footrest (pair)	Std
<input type="checkbox"/>	EL8500Q-16		16" Articulating elevating leg rest with calf support (pair)	\$298.00
<input type="checkbox"/>	EL8500Q-18		18" Articulating elevating leg rest with calf support (pair)	\$298.00
<input type="checkbox"/>	EL8500Q-20		20" Articulating elevating leg rest with calf support (pair)	\$298.00

SEAT CUSHIONS				
✓	SKU#	HCPCS	Description	Price
<input type="checkbox"/>	CU-ERGO-16	E2601	16x17 seat memory foam 2" cushion	\$196.00
<input type="checkbox"/>	CU-ERGO-18	E2601	18x17 seat memory foam 2" cushion	\$176.00
<input type="checkbox"/>	CU-ERGO-20	E2601	20x17 seat memory foam 2" cushion	\$216.00
<input type="checkbox"/>	CU-FT-16	E2601	16x16 seat foam 2" cushion	\$72.00
<input type="checkbox"/>	CU-FT-18	E2601	18x16 seat foam 2" cushion	\$72.00
<input type="checkbox"/>	CU-FT-20	E2601	20x16 seat foam 2" cushion	\$78.00
<input type="checkbox"/>	CU-GFR-16	E2601	16x16 seat gel 2" cushion	\$98.00
<input type="checkbox"/>	CU-GFR-18	E2601	18x16 seat gel cushion	\$98.00
<input type="checkbox"/>	CU-GFR-20	E2601	20x16 seat gel 2" cushion	\$108.00

BACK CUSHIONS				
<input type="checkbox"/>	CU-ERGO-16-back	E2611	16" wide memory foam 2" back cushion	\$156.00
<input type="checkbox"/>	CU-ERGO-18-back	E2611	18" wide memory foam 2" back cushion	\$156.00
<input type="checkbox"/>	CU-ERGO-20-back	E2611	20" wide memory foam 2" back cushion	\$176.00
<input type="checkbox"/>	BKF-1617	E2611	16" wide back foam cushion	\$98.00
<input type="checkbox"/>	BKF-1817	E2611	18" wide back foam cushion	\$98.00
<input type="checkbox"/>	BKF-2019	E2611	20" wide back foam cushion	\$126.00

OPTIONS & ACCESSORIES				
<input type="checkbox"/>	SB22	E0978	Clamp style seat belt 2 piece	\$19.00
<input type="checkbox"/>	SB88	E0978	Wrap around seat belt 1 piece	N/A
<input type="checkbox"/>	SB99-48	E0978	Auto Style seat belt 2 piece	\$24.00
<input type="checkbox"/>	BT-115	E0961	Brake tip-Extending	Std
<input type="checkbox"/>	HR-FLD-115	E0955 + E1028	16-18" width Foldable Headrest with hardware	\$327.00
<input type="checkbox"/>	HR-FLD-115W	E0955 + E1028	20-22" width Foldable Headrest with hardware	\$357.00
<input type="checkbox"/>	AT-115	E0971	Anti Tippers-pair	\$100.00
<input type="checkbox"/>	LS16	K0038	16" width leg strap-calf support	Std
<input type="checkbox"/>	LS18	K0038	18" width leg strap-calf support	Std
<input type="checkbox"/>	LS20	K0038	20" width leg strap-calf support	Std
<input type="checkbox"/>	BKR-EXT-S-105-16		16" width back rest extension	N/A
<input type="checkbox"/>	BKR-EXT-S-105-18		18" width back rest extension	N/A
<input type="checkbox"/>	PSG25F60-BLK	E1015	Frog Legs front fork suspension pair. Patient weight <160 lb	\$299.00
<input type="checkbox"/>	PSG25F61-BLK	E1015	Frog Legs front fork suspension pair. Patient weight >160 lb	\$299.00
<input type="checkbox"/>	CP2	E1399	Universal carry pouch-small	\$11.00
<input type="checkbox"/>	CP4	E1399	Universal carry pouch-large	\$19.00
<input type="checkbox"/>	CANO-115		UV protectant Canopy for 7/8" handle frame	N/A
<input type="checkbox"/>	PB-115-16		Foldable rigidifying push bar for 16" width	\$88.00
<input type="checkbox"/>	PB-115-18		Foldable rigidifying push bar for 16" width	\$88.00
<input type="checkbox"/>	PB-115-20		Foldable rigidifying push bar for 16" width	\$88.00