

KW-100 Luxury Knee Walker

Lightweight 4-wheel Knee Walker w/ Basket



- Steerable unit for easy use and everyday drivability
- Parking brake and slow down allows for safe stopping
- Hub Brake design is the safest in the market
- 8" solid flat free wheels can be used indoors/outdoors
- Ergonomic hand grips help prevent pressure on hands
- Handles collapse for easy carry and storage
- Tool-free handle height adjustment from 31" - 40.5"
- Tool-free seat height adjustment from 17" - 23.5"
- Pad size: 14" x 6"
- 8" flat-free casters for extra stability
- Overall length: 32"; w/basket 37"
- Overall width: 16.5"
- Weighs only 22.5 lbs.
- Weight capacity: 350 lbs.
- Frame Colors: Black, White & Burgundy



Enclosed direct axle contact braking system to stop on a dime and ensure vibrations do not affect your safety and stand the duration of time.



Reinforced woven carbon style ergonomic support for your knee or to sit by rotation for relaxation. Foldable design.



Thickest S-Shape Frame for reaching ultimate stability and style. Height adjustable seat.



KARMAN HEALTHCARE Inc.

19255 San Jose Avenue City of Industry, CA 91748-1418
 Order.1-800-805-2762 Tel: 626-581-2235 Fax:626-581-2335
www.karmanhealthcare.com